

When Diagnostic Security Is At Stake Choose The Best Training: The Simulator.


No book, no video, no CD or DVD can replace interactive real-time training with a realistic investigation scenario.

Only this way the necessary combination of medical know-how with manual skills can be acquired. Only this way all pathologies including rare cases can be comprisingly trained.



Exam 09 - DETA(0) - 2010
max.affon@bini - max.affon@max

SONOFIT
IMAGINE YOUR KNOW-HOW



The simulator behaves practically like an ultrasound scanner and thus permits the selection of stored patients respectively cases/pathologies. Their investigation is performed by section views which are visualized in real-time according to the position of the abdominal, breast and transvaginal probe. Particularly the latter permits training without limitation due to inherent impairment of the patient. Thus, not only training of students but also training of post-graduate physicians is enabled through pathologies which are rarely accessible due to rare incidence.



The system contains an image data base with case/task description and tutorial software to find the correct section view, to do biometry etc.

It is equally suitable for autonomous training as for supporting teaching events, examination of learning progress and quality assurance.

SONOFIT
IMAGINE YOUR KNOW-HOW

SONOFIT GmbH
Am Kirschgarten 5
D-55271 Stackeden-Elshem
Telefon +49 (0)61 30 - 94 14 12
Telefax +49 (0)61 30 - 94 14 45
www.sonofit.com · info@sonofit.com

TROMOZ
TROMOZ