

# TeleConsult 2.0

TeleConsult is the telemedicine-upgrade of the MedCom-InViVo-Software-family.


Based on the MedCom ScanNT-technology it enables its users to perform online and offline medical consultations.




This means by consultation of additional medical experts medical diagnosis can be stated more precisely and certified.

TeleConsult offers the following features:


- Creating and handling of 2D and 3D-Images, including the acquisition of 2D-3D-ultrasound-data and the import of images in DICOM-format.
- Sending and receiving of messages with text and image material over Internet, ISDN and analog phonenumber in offline and online-mode.
- Enrichment of the images with graphical and textual annotations and pictograms.
- In the online-mode interactive communication via usage of a chat window, mouse, exchanged images and a shared screen view is possible.
- Client/Server-architecture of the database. Storage of all images, patients and messages in the database.
- Handling of a second opinion service for customers and partners of a telemedicine center.




consolidated screen view



TeleConsult user interface



Chat: Interactive communication



Example of TeleConsult used in a telemedical center



MedCom Gesellschaft für medizinische Bildverarbeitung mbH  
Rundeturmstraße 6  
D - 64283 Darmstadt, Germany

Phone: +49 (0) 6151 - 95147 - 0  
Fax: +49 (0) 6151 - 95147 - 20

Email: info@medcom-online.de  
URL: http://www.medcom-online.de